

## Laboratory Centrifuge MPW 56

MPW 56 is a small universal centrifuge, dedicated especially for family doctors and veterinaries.

Its structure guarantees a wide range of applications in medical, veterinary, biochemical and other laboratories where constituents of investigated material are to be separated under the effect of centrifugal force.

There are a wide range of angle rotors for tubes with a capacity from 0.2ml to 15ml as well as a haematocrit rotor.



### ADVANTAGES

- Maintenance free brushless motor
- Small size with max. capacity
- Simple control system
- Time pre-selection up to 60 minutes
- Short time of acceleration & deceleration
- Digital LCD display of work time
- Won't start when the lid is open
- Can't open lid during spinning

### Technical Data

Electric supply:	240V 50/60Hz
Power consumption:	70W
Max capacity:	90ml
Speed:	100-6000 rpm
Max RCF:	3340 xg
Time:	0-60 min (15 sec interval adjustment)
Dimensions:	- depth 270mm - width 220mm - height 180mm
Weight without accessories :	5kg

### Other Centrifuges Available From MPW Include:

- MPW 54/55 Small Lab Micro Centrifuges
- MPW 65R Desktop Lab Refrigerated Micro Centrifuge
- MPW 215 Lab Centrifuge
- MPW 223e Table Top Lab centrifuge for Column Tests
- MPW 251 Universal Table Top Lab Centrifuge
- MPW 260 Ventilated / MPW 260R Refrigerated Universal Table Top Centrifuges
- MPW 350, 350R, 350H Universal Table Top Lab Centrifuges
- MPW 351E Table Top Lab Centrifuge

### FOR FURTHER INFORMATION CONTACT:

NSW 02 96091755, VIC 03 9364 9569 or QLD 07 3711 9262

Website: [www.hdscientific.com.au](http://www.hdscientific.com.au)