

Pathology

Centrifuge Special Offer



EBA 20 with rotor
Cat No HE 2002

EBA 20 - Ideal for Collection Centres, Haematology and Blood Bank Depts.

- Complete with 8 place angle rotor for 10ml blood tubes and smaller blood tubes
- Small footprint (around A4 page)
- Weight approx 4 kilos, easy to move around
- Quiet operation, less noise in collection area
- Max speed 6,000 rpm for good separation
- Automatic lid locking conforming to NATA requirements

Regular Price for centrifuge and rotor \$2,086.00

Special Price at \$1,790.00



Rotofix 32A
Cat No HE 1206

4 PI Swing Out Rotor
Cat No HE 1624

4 Buckets
Cat No HE 1381

4 Lids
Cat No HE 1382

4 Inserts
Cat No HE 1329

The Rotofix 32A – A Blood Collection Centrifuge with complete containment.

- Simple operation
- Bio-containment with sealed buckets
- Accepts a wide range of blood tubes up to 10ml.
- Gas strut to protect lid hinges and locking mechanism

Regular Price is including listed accessories \$6,280.00

Special Price at \$5,400.00



Rotolavit
Cat No HE 1006

12 Place Rotor
Cat No HE 1017A

Rotolavit – Automatic Cell Washing System for Serological Testing.

- 12 place or 24 place exchangeable rotors
- Effortless cleaning
- Wash cycle repeated automatically
- Up to 9 sequenced and programmable cycles
- Visual and acoustic signal warns when the saline solution container is empty

Regular Price for centrifuge and 12 place rotor \$21,165.00

Special price at \$18,200.00



Rotina 420
Cat No HE 4701

4 PI Swing Out Rotor
Cat No HE 4723

4 Buckets
Cat No HE 4750

4 Lids with bio containment
Cat No HE 4751

4 Inserts for 10ml Blood tubes
Cat No HE 4763

Rotina 420 – Developed for large volumes, yet it comes with small footprint. High performance when it counts, yet cost-effective and kind to your budget.

- High capacity – 72 x 10ml blood tubes
- Accepts a wide range of other tubes
- Fast run up, 41 seconds to maximum speed
- Whisper quiet operation
- Complete Reliability

Regular Price is \$15,040.00

Special Price at \$13,900.00

Contact us for special deals on Refrigerated Models



Rotina 380
Cat No HE 1701

4 Place swing out rotor
Cat No HE 1754

4 Buckets
Cat No HE 1752

4 Lids with bio containment
Cat No HE 1751

4 Insert 10ml blood tubes
Cat No HE 1763

Rotina 380 – High performance, high capacity and inexpensive.

- Fast “run up” and “run down” times
- Capacity 52 x 10ml blood tubes
- Accepts a wide range of other tubes
- Simple to operate
- Powered lid locking protecting lid lock and hinges
- Programmable

Regular Price is \$13,054.00

Special Price at \$11,753.00

Contact us for special deals on Refrigerated Models

All Hettich Centrifuges come with a 2 year warranty. Preventative Maintenance packages are available for Centrifuges.

For further information complete the coupon below and fax to (02) 9609 1411 or visit our website at www.hdscientific.com.au

Sales Offices: Sydney, Melbourne and Brisbane

Surname:		First Name:	
Title/Position:		Organization Name:	
Address:			
Town/Suburb/State/Postcode:			
Ph:	Fx:	Email:	
Require Information On Service / Preventative Maintenance Packages			
Products Requiring Additional Information			
EBA 20		Rotolavit	
Rotofix 32A		Rotina 420	
Rotina 380			
We are interested in a demonstration of the following model:-			

Conditions apply.

April 2009