



CutMATE

Forceps for Simple Rapid and Consistent Sampling of Specimens at Pre-Set Thickness.

1, 2, 3, 4 An Easy Sequence

1



2



3



4



Advantages of the CutMATE:

- Eliminates the guess work in cutting specific block thickness
- Allows standardization of block thickness, independent of operator ability
- Allows taking serial wet tissue blocks in one simple step
- Suitable for block preparation for use in both microwave and conventional histoprocessors
- Improves confidence and productivity in grossing

For Further Information Contact:

NSW 02 96091755, VIC 03 9364 9569 or QLD 07 3711 9262

Website: www.hdscientific.com.au